



## SUZUKI REACTION MONITORING USING COMPACT MASS SPECTROMETRY WITH THE PLATE EXPRESS™ TLC PLATE READER

### APPLICATION NOTE

**Mass Spec:** [expression](#) CMS

**Sampling:** Plate Express


*Advion, Inc.*

This application note presents an on-line TLC/CMS technique using the Advion [expression](#) Compact Mass Spectrometer (CMS) and Plate Express interface to provide compound structural information without sample preparation after TLC separation for a Suzuki reaction for the synthesis of 4-aminobiphenyl.

## INTRODUCTION

Thin layer chromatography (TLC) is used in many organic synthetic and medical laboratories because it is a simple, cost-effective technique that provides chemists with critical information about their synthetic reactions. Structural characterization of the analytes by TLC is not possible by optical methods such as UV or ELSD. Typically, characterization is performed by GC/MS or LC/MS using sample preparation techniques which involves scraping the TLC spot of interest, extraction using suitable solvents, concentration, and then reconstitution in MS appropriate solvents.


This application note presents an on-line TLC/CMS technique using the Advion *expression* CMS and Advion Plate Express to provide compound structural information without sample preparation after TLC separation. A Suzuki reaction for the synthesis of 4-aminobiphenyl will be demonstrated. The online TLC/CMS technique provides rapid and accurate determination of reaction mixture components without the need for off-line TLC sample preparation procedures.



**Figure 1:** Instrument setup with the CMS and Plate Express.

## METHOD

Reactants A and B were mixed at equimolar amounts in a round-bottom reaction flask and stirred at room temperature. 2 mL aliquots were transferred from the flask and spotted onto a Merck TLC Silica gel 60 F254 plate (10 x 10cm). Chemicals were purchased from Sigma-Aldrich with a purity greater than 99%. MS solvent was LC/MS grade.



**Figure 2:** Chemical structure and information of reactants and catalysts for the Suzuki reaction synthesis of 4-aminobiphenyl.

### TLC Plate Preparation

Merck TLC Silica gel 60 F<sub>254</sub> on Aluminum, 20 x 20 cm - cut to 10 x 10 cm. TLC plates were baked at 100 °C for 10 minutes to displace moisture. A 2  $\mu$ L aliquot of the reaction mixture was spotted onto the TLC plate at the position of 1 cm away from bottom. The TLC plate was then baked at 80 °C for 5 min. 30 mL Benzene was added to a TLC plate developing chamber at the level of 0.5 cm. The developing chamber was sealed with a glass lid for 30 min. The separation was stopped when the development solvent front reached a position which was 1 cm away from the top of the TLC plate.

The analytes (reactant and product) on the developed plate were observed under UV at 254 nm (Figure 3). The R<sub>f</sub> value of the product (4-aminobiphenyl) was 0.1; the R<sub>f</sub> of the reactant (4-bromoaniline) was 0.15.




Figure 3: The developed TLC plate with Suzuki reaction mixture.

### TLC/CMS Method

The TLC/CMS analysis of the Suzuki reaction mixture at different reaction times was performed on the Advion **expression** CMS and Plate Express. A solvent composed of 0.1% formic in methanol was used for the elution of the analytes from the TLC plate. The eluted analytes were directed to the CMS for acquisition of the corresponding mass spectra for the reactants and products.

## RESULTS

The mass spectra of the reactant and product are shown in Figure 4. The [M+H]<sup>+</sup> for 4-bromoaniline was observed at  $m/z$  171.9 (<sup>79</sup>Br) and 173.9 (<sup>81</sup>Br). Methanol adducts were also detected at  $m/z$  204.0 and 206.0 (Figure 4A).

The [M+H]<sup>+</sup> for 4-aminobiphenyl was observed at  $m/z$  170.1 with a methanol adduct at  $m/z$  202.1 (Notice the absence of the bromine doublet peaks in the non-brominated product of p-aminobiphenyl, Figure 2B).





Figure 4: TLC/CMS analysis of reactant and product. (A) The mass spectrum of protonated 4-bromoaniline (reactant) at the  $m/z$  171.9 and 173.9. (B) The mass spectrum of protonated 4-aminobiphenyl (product) at the  $m/z$  170.1

The extracted ion current of the product ion at  $m/z$  170.1 (protonated 4-aminobiphenyl) was monitored over different reaction times from 0 to 180 minutes (Figure 5A). The mass spectrum of the product ion at 120 min is shown in Figure 5B. The extracted ion current (XIC) of the reactant ion at the  $m/z$  171.9 (protonated 4-bromoaniline( $^{79}\text{Br}$ )) is shown in Figure 5C. The intensity of the  $m/z$  171.9 decreased over the course of the reaction with respect to the product ion at the  $m/z$  170.1. The reactant ion at the  $m/z$  171.9 was no longer detected at 180 minutes, indicating the reaction was complete.

## SUMMARY

- The Advion **expression** CMS coupled with the Advion Plate Express offers a simple and fast technique to monitor a Suzuki reaction for the synthesis of 4-aminobiphenyl.
- The Advion TLC/CMS system allows the synthetic chemist to monitor the reaction in real-time by evaluating the mass spectra for structural information (i.e., relative intensity of reactants vs. product) directly from the TLC plate.
- The compact size allows it to fit into space-limited labs for direct access and immediate results for chemists requiring mass confirmation, reaction monitoring, quality control and purity analysis.



**Figure 5:** (A) XIC of product ion (protonated 4-aminobiphenyl), (B) mass spectrum of the product ion at reaction time at 120 min, (C) XIC of the reactant ion (protonated 4-bromoaniline), and (D) mass spectrum of the reactant ions at 15 minutes.